

THE ORASTRETCH® PRESS AND THERABITE® DEVICE JAW MOTION REHAB SYSTEMS

The **OraStretch** Press and **TheraBite** jaw motion rehab systems are handheld passive motion devices designed for the prevention and treatment of trismus, facial scarring, surgical rehab, pain and joint dysfunction.



Use of a jaw motion rehab system has been shown to reduce the incidence by and the severity of trismus in cancer patients^{1,2}.

The **OraStretch** press and the **TheraBite** move along the natural rotation of the jaw for easy, safe rehabilitation with limits for control and safety. Both include a case, manual, scales, and pads.



Pediatric versions are available for the **OraStretch** press and **TheraBite** devices for patients under 10 years-old and small mouths. Edentulous pads for missing teeth or in place of dentures are also available for each system.

The **OraStretch** press also offers an extended version for improved dental stability in patients with a severe under- or overjet; and further includes free replacement pads, free size exchanges, and a lifetime warranty.



Contact us for more information on insurance coverage, pricing, ordering, usage, and recommended protocols.

Trismus * Radiation Fibrosis * TMJD * Stroke * Post-Surgical Rehab

Jaw Motion Rehab Systems	OraStretch Press	TheraBite Device
		
Versions	<ul style="list-style-type: none"> • Standard • Pediatric • Extended 	<ul style="list-style-type: none"> • Standard • Pediatric
Therapeutic Motion	5" radius to 50mm	5" radius to 41mm
Weight	3.8 oz	7.7 oz
Mouthpiece Width	Std/Ext: 48mm, Ped: 31mm	Std: 48mm; Ped: 31mm
Safety Limits	Yes: Limit Screw	Yes: Limit Screw & Arm Pin
Scale	Yes: Visible on Side	Yes: Under chin
Warranty	Lifetime	NONE
Advantages	<ul style="list-style-type: none"> • Billable to Medicare. • 20% wider ROM to full adult average. <ul style="list-style-type: none"> • Maximum trismus prevention. • 50% smaller, easier to hold & use. • Visible side scale. • Free Size Exchanges*. • Lifetime Warranty*. 	<ul style="list-style-type: none"> • Recognized brand. • Hand slide for severe arthritis.

1: SONG Dan,JIANG Ming-hua. Effect of man-made intervention on IMRT-induced trismus in patients with nasopharyngeal carcinoma. Chinese Clinical Oncology;2008-11, 2: Tang Y, Shen Q, et al. A randomized prospective study of rehab therapy in the treatment of radiation-induced dysphagia and trismus. Strahlenther Onkol. 2011 Jan;187(1):39-44. 3: Melchers LJ, Van Weert E, et al. Exercise adherence in patients with trismus due to head and neck oncology. Int J Oral Maxillofac Surg. 2009 Sep;38(9):947-54.

* See Website for Details. US Only. Prices and Offers Subject to Change.



CRANIOMANDIBULAR REHAB, INC.
"HELPING PEOPLE SAY AHHHH!"™

WWW.CRANIOREHAB.COM
1-800-206-8381